



PLURYAL INJECTABLES

products range
overview

FILLERS ESSENTIALS



GENERAL PARAMETERS	PLURYAL® BOOSTER	PLURYAL® CLASSIC	PLURYAL® VOLUME
COMPOSITION	HA : 18mg/ml Glycerol : 20mg/ml	Medium Cross-linked HA: 23mg/ml	Highly Cross-linked HA: 23mg/ml
INDICATIONS	Skin elasticity and hydration	Loss of volume Wrinkles	Loss of volume Deep wrinkles
TARGET ZONES	Face Neck Décolletage Back of the hands	Lips Nasolabial folds Tear trough Neck	Nose Deep wrinkles Nasolabial folds Scars Facial lipoatrophy Morphological asymmetry
TECHNOLOGY	PREMIUM	PREMIUM	PREMIUM
HA MOLECULAR WEIGHT (Da)	3 Million	3 Million	3 Million
RETICULATION	-	●●●○○	●●●●○
ELASTIC MODULUS G' (mPa)	-	195 000	400 000
RESIDUAL BDDE (ppm)	-	< 2	< 2
LEVEL OF INJECTION	Superficial dermis	Mid to deep dermis / Submucosa	Deep dermis or subcutis
PACKAGING	1 x 1ml	1 x 1ml	1 x 1ml
CE MARKING	CE0123	CE0123	CE0123

FILLERS LIDOCAINE



GENERAL PARAMETERS	PLURYAL® CLASSIC LIDOCAINE	PLURYAL® VOLUME LIDOCAINE	PLURYAL® CONTOUR LIDOCAINE
COMPOSITION	Medium Cross-linked HA: 23mg/ml with lidocaine 3mg/ml	Highly Cross-linked HA: 23mg/ml with lidocaine 3mg/ml	Very highly Cross-linked HA: 25mg/ml with lidocaine 3mg/ml
INDICATIONS	Loss of volume Wrinkles	Loss of volume Deep wrinkles	Loss of volume and asymmetries
TARGET ZONES	Lips Nasolabial folds Tear trough Neck	Nose Deep wrinkles Nasolabial folds Scars Facial lipoatrophy Morphological asymmetry	Temples Midface volume deficit Jawlines Chin
TECHNOLOGY	PREMIUM	PREMIUM	PREMIUM
HA MOLECULAR WEIGHT (Da)	3 Million	3 Million	3 Million
RETICULATION	●●●○○	●●●●○	●●●●●
ELASTIC MODULUS G' (mPa)	195 000	400 000	650 000
RESIDUAL BDDE (ppm)	< 2	< 2	< 2
LEVEL OF INJECTION	Mid to deep dermis / Submucosa	Deep dermis or subcutis	Deep subcutaneous and/or supraperiosteal
PACKAGING	1 x 1ml	1 x 1ml	1 x 1ml
CE MARKING	CE 0123	CE 0123	CE 0123

DENSIFIER/FILLERS BIO



GENERAL PARAMETERS	PLURYAL® BIOSCULPTURE	PLURYAL® BIOCLASSIC	PLURYAL® BIOVOLUME
COMPOSITION	Low Cross-linked HA & Free HA: 25mg/ml	Medium Cross-linked HA & Free HA: 25mg/ml	Highly Cross-linked HA & Free HA: 25mg/ml
INDICATIONS	Sagging skin Laxity Flaccidity	Loss of volume Wrinkles Mid laxity and dehydration	Loss of volume Wrinkles Deep laxity and dehydration
TARGET ZONES	Face Back of the hands	Lips Nasolabial folds Tear trough Neck Cheeks	Temples Cheeks Chin - Jawline Lips Nasolabial folds
TECHNOLOGY	MATRICYAL 3D	MATRICYAL 3D	MATRICYAL 3D
HA MOLECULAR WEIGHT (Da)	Cross-linked Matrix : 2 Million (60%), 1 Million (30%) Free HA (10%): 0,5 Million	Cross-linked Matrix: 2 Million (78%), 1 Million (15%) Free HA (7%): 1 Million	Cross-linked Matrix: 2 Million (90%), 1 Million (5%) Free HA (5%): 1 Million
RETICULATION	●●○○○	●●●○○	●●●●○
ELASTIC MODULUS G' (mPa)	-	100 000	200 000
RESIDUAL BDDE (ppm)	< 0.1	< 0.1	< 0.1
LEVEL OF INJECTION	Mid to deep dermis	Mid to deep dermis	Hypodermis
PACKAGING	2 x 1,1 ml	1 x 1ml	1 x 1ml
CE MARKING	CE0373	CE0373	CE0373

DENSIFIERS HPN



GENERAL PARAMETERS	PLURYAL® SILK	PLURYAL® DENSIFY	PLURYAL® HAIR DENSITY
COMPOSITION	Polynucleotides : 7,5mg/ml	Polynucleotides : 10 mg/ml HA : 10 mg/ml Mannitol : 40 mg/ml	Polynucleotides : 7,5mg/ml
INDICATIONS	Aggressed skin Skin repair : texture and color Preparation for other aesthetic treatments	Damaged, dehydrated skin Skin repair : texture and color Preparation for other aesthetic treatments	Hair growth Hair density
TARGET ZONES	Delicate areas (eye and mouth contour) Stretch marks, scars	Full face (except around the eyes) Back of the hands Neck Décolletage	Scalp Eyebrows Beard
TECHNOLOGY	HPN	HPN	HPN
PN CHARACTERISTICS	<p>Biopolymers composed of nucleotide monomers</p> <p>Extracted from salmon milt and obtained from a fraction of DNA of different chain lengths (molecular weight from 80 to 2,000 kDA) close to the human DNA</p> <p>Highly purified to remove any traces of protein</p> <p>Non-allergenic and fully absorbable</p> <p>Our DNA is formed from this kind of molecules, making them very important molecules for our body's functioning</p> <p>They have the ability to repair the skin and stimulate the hair growth</p>		
LEVEL OF INJECTION	Deep dermis	Deep dermis, superficial hypodermis	Deep dermis
PACKAGING	1 x 2ml	1 x 2ml	1 x 2ml
CE MARKING	CE0373	CE0373	CE0373

YOUR MOST BEAUTIFUL AGE IS NOW!



Pluryal is a registered trademark of MD Skin Solutions.

Pluryal
MD Skin Solutions
9AB Boulevard du Prince Henri
L-1724 Luxembourg
T: +352 26 20 04 06
pluryal.com

P_LTES_PRO_EN_RevE